# Instrument Cleaning Care CLARINET

Online Video at: www.rockymountainmusicrepair.com/instrument-care

### Daily

- Remove water from inside of the clarinet with a drop swab.
   Remove mouthpiece, put weight in from bell, tilt clarinet so that the weight comes out the barrel and pull though. (If it ever gets stuck, take it to a repair shop right away, don't try pulling on it harder).
- Store reeds in reed holder or case, NEVER on mouthpiece.

### Weekly

- Gently remove fingerprints from the keys with soft cloth.
- Clean mouthpiece with a mouthpiece brush and warm soapy water.
- Clean & grease all tenon corks with a dab of new cork grease.
- Vacuum case to get rid of all debris.

### Yearly

• Bring your instrument to a qualified repair technician for a professional regulation & lubrication and to replace needed pads.

## **Cleaning Supplies Needed**

 Clarinet Drop Swab, Soft Cloth, Woodwind Mouthpiece Brush, Cork Grease



